

Product PORTFOLIO

GRIFOLS

At Grifols, we make it a priority to listen and learn about our customers and the patients they care for. This understanding informs all that we do, and fosters the partnerships that enrich our work and our goals for the future. Grifols USA Hospital Division provides simple, cost-effective solutions to operational challenges in pharmacy. Its primary focus is on sterile IV preparations for specialty therapies such as pain management, oncology, and pediatrics.



Gri-fill® 3.0 – Gri-fill® is a closed system which allows sterile solutions to be prepared in different volumes, contents, and proportions than those produced by the pharmaceutical industry, all without the need for laminar flow hoods. The Gri-fill® 3.0 System can be placed inside a standard laminar air flow hood or biosafety cabinet when extra protection is needed (i.e. with antineoplastics). The Gri-fill® maintains sterility via the principle of sterile filtration and quality control when Gri-bags® and Gri-flex®


containers are used. Each unit undergoes final filtration through a 0.22 micron filter and its integrity is automatically verified. The Gri-fill® System, in accordance with USP Chapter 797, is designed to achieve compounding accuracy and sterility to ensure the safety of patients.



MixManager – MixManager is a Windows®-based software that is part of the Gri-fill® System that organizes your workflow and simplifies the order entry process while allowing you to link multiple compounders for maximum efficiency. Customization of reports and labels help meet documentation requirements and streamline the workflow in compounding areas.



Gri-bags® – The Gri-fill® System bags have the capability to filter every dose, offering increased confidence in delivery of a sterile preparation to your patients. The integral 0.22 micron filter on Gri-bag® containers assures the maintenance of sterility in a closed system backed by an integrity test. The Gri-bag® AP, with Additive Port, offers the flexibility to use multiple set spike sizes or to introduce another additive just prior to administration.



Flebaset® – Streamline your workflow and improve efficiency with a standard set up and standard process. One set is used for each drug-diluent combination. The same set can be used for preparing both batch and custom doses using the same components. The multi-spike set for use with multiple vials is used to create a virtual pooled container, combining multiple glass vials of the same drug concentration. It gives you the advantage of

reducing handling, minimizing the risk of contamination, accommodates multiple vial sizes, and allows vial replacement as needed without interruption of compounding.



MISTERIUM – MINI>STERILE>MODULAR>ROOM...COMING SOON

A modular clean room configured to suit your needs, installed wherever you wish, without the need for expensive building projects, without overshooting budgets or delays in delivery, thanks to its easy-assembly, modular panel system.

Call toll free at 1-800-379-0957 or visit the website at www.grifolsusa.com