

Disinfectant-Ready Mobile Computer

Honeywell announces its Dolphin 9900hc mobile computer with disinfectant-ready housing, specifically designed to mobilize patient safety in environments where continuous cleaning of the



device is required. A mobile computing solution for point-of-care applications, the Dolphin 9900hc's housing resists the harmful effects of repeated exposure to harsh cleaning solvents. The Dolphin 9900hc retains all the features and functions of the Dolphin 9900, including Adaptus Imaging Technology. The computing architecture is powered by Microsoft Windows Mobile 6 and delivers wireless communication that supports advanced wireless security standards.

- **Honeywell Scanning & Mobility**
Circle reader service number 1
 or visit www.honeywell.com/aidc

Oral Solid Feeder

Medical Packaging Inc (MPI) introduces its Auto-Print Oral Solid Feeder featuring FDA-approved materials and coatings. This rotary feeder adapts



to MPI's Auto-Print, streamlining the packaging process for oral solids into unit dose. Bulk medication is loaded into the feeder packaging one dose at a time. When medication supply is depleted, the Auto-Print will automatically pause and notify the user with an audible alert. The Auto-Print Oral Solid Feeder will dispense tablets up to 1/2 inch diameter, capsules, and a variety of other medications.

- **Medical Packaging Inc**
Circle reader service number 2
 or visit www.medpak.com

Pharmacy Pricing Software

Amerinet announces its AccuPrice software tool that combines master file maintenance, daily price audits, and staff resources to identify and resolve wholesaler pricing errors. Each day, the AccuPrice software reviews all purchases and submits an audit report by e-mail to the facility including elements such as unfavorable and



favorable variances on base- as well as tier-pricing programs. Additionally, facility summary reports that capture and detail all pricing errors for specified time frames are compiled and distributed. Because errors are identified nearly immediately, billing staff can contact the supplier or distributor before the errors repeat, thereby decreasing the risk of recurrence. The monthly AccuPrice summary report also provides a single snapshot to track variances along with credits and rebills.

- **Amerinet**
Circle reader service number 3
 or visit www.amerinet-gpo.com

Anticoagulation Alerts Package

TheraDoc introduces a patient safety surveillance tool to help hospitals address the recent Joint Commission National Patient Safety Goal focused on reducing the possibility of patient harm associated with anticoagulation therapy. The TheraDoc Anticoagulation Alerts package

Alert	Alert	Alert	Alert	Alert	Alert
Heparin, new aPTT	Admit Diagnosis: CAD RECENT HF	This patient has one or more active orders for heparin and has a new Partial Thromboplastin Time result available for review.			
Order/Culture	Result	Source	Collected	Result Status (Date/Time)	Specimen #
Partial Thromboplastin Time (Partial Thromboplastin Time = 100 sec; n=1 (24-25)) [See ESP (05/03/2009 20:29) # (05/03/2009 21:31) (04097166)]					
Active or Hold Heparin Orders:					
Drug	Dose	Start	End	Status	Pat Class
HEPARIN	25000 unit IV	05/02/2009 02:31:00		ACTIVE	I
Ten Most Recent Coagulation Modifiers Orders:					
Drug	Dose	Start	End	Status	Pat Class
ASPIRIN	162 mg PO QDay	04/28/2009 08:00:00		ACTIVE	I
CLONIDINE	75 mg PO QDay	04/28/2009 08:00:00	05/03/2009 08:47:00	INACTIVE	I
EPITIFIBATIDE	75 mg IV	05/03/2009 12:00:00		ACTIVE	I
HEPARIN	25000 unit IV TID	04/28/2009 08:00:00	04/28/2009 15:33:44	INACTIVE	I
HEPARIN	25000 unit IV	04/28/2009 15:33:00	05/03/2009 08:48:47	INACTIVE	I
HEPARIN	25000 unit IV	05/03/2009 02:31:00		ACTIVE	I

tracks the condition of patients receiving anticoagulation drug therapy to prevent blood clots, providing an early warning of potential adverse drug events or other problems. Designed for use with the TheraDoc Expert System Platform, the package includes approximately 43 alerts from existing TheraDoc products, and 21 new alerts are scheduled for release in the Anticoagulation Alerts package at the end of June.

- **TheraDoc**
Circle reader service number 4
 or visit www.theradoc.com

Disinfectant

Haddad-Wylie Industries (HWI Global) introduces its Bio-Gard Neutral Quat 64 disinfectant. Bio-Gard is effective against a broad spectrum of bacteria (MRSA, Staphylococcus), is virucidal (Avian, HEP, HIV, Influenza), fungicidal (Candida, Aspergillus), and inhibits the growth and odors of mold and mildew when used as directed. Bio-Gard cleans, disinfects, and deodorizes hard, non-porous environmental medical surfaces in one step with no rinsing required. Starter kit includes 16 oz of Neutral-Quat 64 concentrate and a 32 oz refillable spray bottle. The concentrate will yield 32, 32 oz bottles of ready-to-use Bio-Gard Neutral Quat 64.



- **Haddad-Wylie Industries**
Circle reader service number 5
 or visit www.hwicleanrooms.com

Tablet-Splitting Device

TCGRx (TCG) announces its BullsEye Tablet Splitter, a device that splits medication with two percent variability between the split halves. The



BullsEye Tablet Splitter can automatically split up to 16 pills per minute and contains a replaceable blade and HEPA filter to decrease potential contamination.

- **TCGRx**
Circle reader service number 6
 or visit www.tcgrx.com

Bar Code Scanner

Opticon announces its M5 2D imager. The bar code scanner features auto sensing optics with three different scan modes customized for user functionality. Choose between the auto-trigger mode for hands-free scanning, manual mode for your specific application, or mobile mode, which allows bar codes to be read from the screen of a mobile device. The M5 comes with a built-in adjustable stand allowing adaptation to specific working environments.



► **Opticon Inc**
Circle reader service number 7
 or visit www.opticonusa.com

Ready to Use IV Product

B. Braun announces the DUPLEX drug delivery system. Consisting of a pre-filled, DEHP-free, PVC-free, latex-free IV bag, containing proper doses of drug and diluent in different compartments separated by a quick-release seal, DUPLEX can be stored at room temperature and does not require thawing. Oxygen and moisture barriers are employed to protect the drug



powder from vapor and oxygen transfer through the bag over the product's shelf life. Equipped with a HIBC standard bar code that references the final admixture, the DUPLEX system automates patient charting, tracks inventory, and facilitates reimbursement tracking.

► **B. Braun Medical Inc**
Circle reader service number 8
 or visit www.bbraunusa.com



*Leading
the Way...*

in contamination control for
USP <797> compliance

USP <797> describes the cleaning practice for compounding sterile preparations. ITW Texwipe® provides the cleaning products for compliance:

- Isolator cleaning systems
- Sterile 70% IPA
- Wipers and swabs
- Mopping systems

USP <797> requires cleaning SOPs, and ITW Texwipe offers cleaning guidelines for:

- Isolator cleaning
- Wiping
- Mopping

Protect your investment with the contamination control leadership of ITW Texwipe.

**Sterile IPA • Dry and IPA-wetted Sterile Wipers
 Cleanroom Wipers • Mops • Stationery • Swabs**



	North America	Europe	Asia
Tel	201 684 1800	+45 87 400 220	+65 6468 9433
Fax	201 684 1801	+45 87 400 222	+65 6468 6772
E-mail	info@texwipe.com	europa@texwipe.com	asia@texwipe.com

ITW Texwipe®

Quality. Consistency. Support.

www.texwipe.com

For more information, circle #11 on the Reader Service Card

New

& NOTEWORTHY

Plasma and Biologics Unit

McKesson announces the launch of its Plasma and BioLogics unit, which provides distribution service for plasma and biologics products for pharmaceuticals directly from the manufacturer to customers across the continuum of care. Services include monthly reporting and bill pay-

MCKESSON

Empowering Healthcare

ment, McKesson ordering-portal availability, contribution towards total purchasing volume, dedicated customer service and sales representatives and an emergency order system operating outside regular business hours through a toll-free number.

- ▶ **McKesson Corporation**
Circle reader service number 9
or visit www.mckesson.com

Wireless Temperature Monitoring

Comark's new RF500 system provides a full audit trail for 24/7 monitoring of temperature and humidity. A network of remote sensors and probes collects and transmits information for data management and storage. In the event of a power failure, the transmitters store and send data as soon as the network is restored. A mesh network enables the system to automatically reroute signals as required to maintain integrity of the data.



- ▶ **Comark Instruments**
Circle reader service number 10
or visit www.comarkusa.com

Counting, Verification, and Storage

Kirby Lester announces an alliance with Uniweb to provide an all-in-one automated dispensing



and verification system. Uniweb has developed shelving units that accommodate Kirby Lester machines and accessories. The Uniweb KL30 storage system provides a vertical method to organize 30 or more, drug-specific dispensing cassettes. The wall-mount shelf holds the Kirby Lester tablet counter on an open wall area.

- ▶ **Kirby Lester, LLC**
Circle reader service number 12
or visit www.kirbylester.com

Controlled Environment Mops

Perfex introduces three new mopping products to its product line. The TruCLEAN clean-room



mop (#22-39) is constructed of continuously woven polyester fabric, individually laundered and packaged. The TruCLEAN micro-fiber mop (#22-31) contains many strands of a synthetic fiber made from nylon and a polyester/polyamide blend. The TruCLEAN anti-microbial micro-fiber mop has an anti-microbial agent permanently bonded to the micro-fibers. TruCLEAN mops are ideal for cleaning and disinfecting ceilings, walls, and floors. These mops are compatible with all TruCLEAN mopping systems.

- ▶ **Perfex, Inc**
Circle reader service number 13
or visit www.perfexonline.com

High-Speed Packager

Through an agreement with TCGRx, Swisslog announces the ATP-71, a compact and affordable high-speed packager. The ATP-71 can accommodate up to 71 canisters in a small footprint while offering 100% bar coding for oral solids through a universal feeder. The packager can handle both unit dose bar codes and box replenishment with a packaging speed of up to 60 packages per minute. An external tray is available for the system's non-formulary medications.



- ▶ **Swisslog**
Circle reader service number 14
or visit www.swisslog.com

Advertisers Index:

COMPANY	PAGE	READER SERVICE #	COMPANY	PAGE	READER SERVICE #
The Baker Company	Inside Front Cover	45	NuAire, Inc	13	68
Cooper-Atkins	Back Cover	49	Pevco	9	95
ITW Texwipe	23	11	R.C. Smith Company	Inside Back Cover	72
McKesson Corporation	19	59	Rees Scientific Corporation	8	55
medDISPENSE	21	39	Talyst, Inc	7	32
Medi-Dose, Inc/EPS, Inc	3	15	Terra Universal, Inc	11	77
Medical Packaging, Inc	17	65	Teva Pharmaceuticals	5	100
Modular Cleanrooms Inc	15	85			